

chronic_ACH_exposition_metrics

June 19, 2025

0.1 Control Dataset Overview

This notebook analyzes spontaneous calcium activity in BxPC3 pancreatic cancer cells imaged under identical control conditions across six image sequences.

Experimental context: - Cell line: BxPC3 (pancreatic ductal adenocarcinoma) - Imaging: DIV2, seeded at 1.2 million cells per well (high confluence) - Ca^{2+} dye: Fluo-4 AM (standard loading protocol) - Nuclei: Hoechst 33342 (10 $\mu\text{g/mL}$) - Acquisition: FITC, 1 Hz, 30 ms exposure for 30 minutes - CHRONIC exposition to Ach in growth medium after seeding, during incubation & imaging at various concentration.

This is the **first run** of spontaneous activity under control conditions. The goal is to analyze cells' behaviors at different level and compare it to control datasets

0.2 Setup & Utilities

Standard scientific libraries (Pandas, NumPy, Seaborn) and custom plotting utilities are loaded. These functions enable streamlined metric visualization and layout control across multiple datasets.

```
[1]: # Libraries
import seaborn as sns
import numpy as np
from pathlib import Path
import sys

# Add /Source to sys.path if not already present
project_root = Path.cwd().parent
source_path = project_root / "src"

if source_path.exists() and str(source_path) not in sys.path:
    sys.path.insert(0, str(source_path))

# Utility functions
from calcium_activity_characterization.analysis.visualizers import
    ↪plot_metric_by_dataset, plot_category_distribution_by_dataset,
    ↪plot_raster_plots_by_dataset
from calcium_activity_characterization.analysis.loaders import
    ↪load_and_merge_datasets

sns.set(style="whitegrid")
```

```

# Define paths to control datasets
control_paths = {
    "1 - ACH 10uM - 1": "D:/Mateo/20250618/Output/IS1",
    "1 - ACH 10uM - 2": "D:/Mateo/20250618/Output/IS2",
    "2 - ACH 1uM - 1": "D:/Mateo/20250618/Output/IS3",
    "2 - ACH 1uM - 2": "D:/Mateo/20250618/Output/IS4",
    "3 - ACH 100nM - 1": "D:/Mateo/20250618/Output/IS5",
    "3 - ACH 100nM - 2": "D:/Mateo/20250618/Output/IS6",
    "4 - ACH 10nM - 1": "D:/Mateo/20250618/Output/IS7",
    "4 - ACH 10nM - 2": "D:/Mateo/20250618/Output/IS8",
}

# Load and merge datasets
peaks = load_and_merge_datasets(control_paths=control_paths, file_name="peaks.
↪csv")
cells = load_and_merge_datasets(control_paths=control_paths, file_name="cells.
↪csv")
events = load_and_merge_datasets(control_paths=control_paths, file_name="events.
↪csv")

```

0.3 Preprocessing

Before visual comparison, we normalize event labels (e.g., missing `in_event` entries set to 'individual') and compute additional derived fields (`is_active`, `periodicity_score`) to ensure consistent comparisons across all datasets.

These steps are crucial for **ensuring fair metric aggregation** and enabling downstream comparative analysis.

```

[2]: peaks["in_event"] = peaks["in_event"].fillna("individual").str.lower()

if "is_active" not in cells.columns:
    cells["is_active"] = cells["num_peaks"].apply(lambda x: 1 if x > 0 else 0)

if "periodicity_score" not in cells.columns:
    cells["periodicity_score"] = np.nan

```

0.4 Raster Plot Inspection

Raster plots provide a binary overview of calcium activity over time per cell.

Remarks: TODO

Use these plots to **screen visually for anomalies or batch effects** before statistical analysis.


```

[3]: plot_raster_plots_by_dataset(
    dataset_paths=control_paths,
    title="Binary Activity Raster Plots by Dataset",

```

```
n_cols=2  
)
```

Binary Activity Raster Plots by Dataset



0.5 Peak Type Distribution

Pie charts quantify the proportion of detected calcium peaks classified as: - **individual**: isolated peaks - **sequential**: propagating local events - **global**: large-scale synchronized events

Interpretation:

The distribution of peaks by event type (global, sequential, individual) across First Run control datasets is highly consistent.

- *Approximately **30-40% global events** are consistently observed, except in image sequence 2025-04-04 IS3, where fewer global events were detected.*
- *Sequential and individual peaks maintain stable proportions across datasets, indicating reproducible spontaneous calcium signaling activity under standard conditions.*

First run image sequences demonstrate reproducible and robust peak-type distribution, validating their reliability as control data.

```
[4]: peak_type_colors = {  
    "global": "#66c2a5",  
    "sequential": "#fc8d62",  
    "individual": "#8da0cb"  
}  
  
plot_category_distribution_by_dataset(  
    df=peaks,  
    column="in_event",  
    category_order=["global", "sequential", "individual"],  
    colors=peak_type_colors,  
    title="Distribution of Peaks by Event Type",  
    n_cols=2  
)
```

Distribution of Peaks by Event Type



0.6 Peak Metrics per Event Type

Metrics like **duration**, **prominence**, and **symmetry** are visualized per event type.

Interpretation:


- *Global peaks: Characteristically long durations (~50 frames).*
- *Sequential peaks: Duration means around **26–28 frames**; symmetry values around **0.9 ± 0.12**.*
- *Individual peaks: Duration around **26–27 frames**, also with symmetry $\sim 0.9 \pm 0.12$.*

Overall, the peak metrics in the First Run are stable, biologically plausible, and reproducible.


```
[5]: for event_type in ["global", "sequential", "individual"]:
      subset = peaks[peaks["in_event"] == event_type]
      print(f"\n# Peak Type: {event_type.title()} - {len(subset)} peaks")
      plot_metric_by_dataset(subset, "rel_duration", f"{event_type.title()} Peaks:
↪ Duration", bin_width=3, n_cols=2)
      plot_metric_by_dataset(subset, "prominence", f"{event_type.title()} Peaks:
↪ Prominence", bin_width=3, n_cols=2)
      plot_metric_by_dataset(subset, "rel_symmetry_score", f"{event_type.title()}
↪ Peaks: Symmetry", bin_count=30, n_cols=2)
```

```
# Peak Type: Global - 35737 peaks
```


Global Peaks: Duration



Global Peaks: Prominence




Global Peaks: Symmetry




Peak Type: Sequential - 17127 peaks


Sequential Peaks: Duration



Sequential Peaks: Prominence




Sequential Peaks: Symmetry




Peak Type: Individual - 9321 peaks


Individual Peaks: Duration



Individual Peaks: Prominence



Individual Peaks: Symmetry



0.7 Cell Metrics

We analyze: - Number of peaks per cell (activity) - Periodicity score (oscillatory behavior) - Proportion of active vs inactive cells

Interpretation:

Active cells: Nearly all cells (>99.9%) are active in every dataset. This indicates excellent signal integrity and robust experimental conditions. Activity is highly homogeneous across datasets.









Number of peaks per cell: High total peak count (>10,000 peaks per dataset), with a balanced and near-normal distribution. This reflects widespread and stable activity across the population.

Periodicity: Periodicity scores ($\sim 0.7 \pm 0.1$) remain consistent across both imaging runs. This suggests that rhythmic calcium oscillations are preserved, even when global excitability drops.

```
[6]: # Convert is_active to label for display purposes
cells["active_label"] = cells["is_active"].map({1: "Active", 0: "Inactive"})


plot_category_distribution_by_dataset(
    df=cells,
    column="active_label",
    category_order=["Active", "Inactive"],
    colors={"Active": "#66c2a5", "Inactive": "#d3d3d3"},
    title="Active vs Inactive Cells per Dataset",
    n_cols=2
)
```

Active vs Inactive Cells per Dataset




```
[7]: plot_metric_by_dataset(cells, "num_peaks", "Cell: Number of Peaks",  
    ↪bin_width=1, n_cols=2)  
plot_metric_by_dataset(cells, "periodicity_score", "Cell: Periodicity Score",  
    ↪bin_width=0.02, n_cols=2)
```

Cell: Number of Peaks



Cell: Periodicity Score



0.8 Sequential Event Metrics

We focus on **spatiotemporal events**, extracting: - Mean communication time and speed - Geometric descriptors (elongation, radiality) - Graph depth, number of involved cells


Interpretation:

- *Communication time: ~ 4.0 frames (STD ~ 2.5).*
- *Communication speed: ~ 6.2 units/frame.*
- *Cells involved: $\sim 3.2 \pm 2$ cells per event.*


These metrics are consistent across datasets and define the propagation baseline for healthy spontaneous events.

```
[8]: seq_events = events[events["event_type"] == "SequentialEvent"]
plot_metric_by_dataset(seq_events, "communication_time_mean", "Sequential_
↳Events: Comm Time", bin_width=1, n_cols=2)
plot_metric_by_dataset(seq_events, "communication_speed_mean", "Sequential_
↳Events: Comm Speed", bin_count=50, n_cols=2)
plot_metric_by_dataset(seq_events, "elongation_score", "Sequential Events:
↳Elongation", bin_count=50, n_cols=2)
plot_metric_by_dataset(seq_events, "radiality_score", "Sequential Events:
↳Radiality", bin_width=0.05, n_cols=2)
plot_metric_by_dataset(seq_events, "dag_depth", "Sequential Events: DAG Depth",
↳bin_width=1, n_cols=2)
plot_metric_by_dataset(seq_events, "n_cells_involved", "Sequential Events:
↳Cells Involved", bin_width=1, n_cols=2)
```


Sequential Events: Comm Time




Sequential Events: Comm Speed



Sequential Events: Elongation



Sequential Events: Radiality



Sequential Events: DAG Depth

1 - ACH 10uM - 1 (No Data)

1 - ACH 10uM - 2 (No Data)

2 - ACH 1uM - 1 (No Data)

2 - ACH 1uM - 2 (No Data)


3 - ACH 100nM - 1 (No Data)

3 - ACH 100nM - 2 (No Data)

4 - ACH 10nM - 1 (No Data)

4 - ACH 10nM - 2 (No Data)

Sequential Events: Cells Involved



0.9 Conclusion & Next Steps

0.9.1 Summary of Observations

- **Visual Patterns:** Briefly describe observed visual trends (raster plots, activity levels, etc.).
- **Statistical Results:** Mention key metrics or KS-test outcomes that stood out.

0.9.2 Dataset Reproducibility

- Which datasets appear consistent or divergent?
- Are any replicates potentially faulty or biologically distinct?

0.9.3 Actionable Decisions

- ☐ Exclude dataset(s) due to technical issues
- ☐ Flag dataset(s) for further inspection
- ☐ Select representative samples for downstream analysis

0.9.4 Planned Deep-Dive

- Which direction will the next notebook take? (e.g., spatial clustering, propagation analysis, signal decomposition)

```
[ ]: from calcium_activity_characterization.analysis.report import ↪export_current_notebook_to_pdf
#export_current_notebook_to_pdf("chronic_ACH_exposition_metrics.ipynb")
```

Notebook export error: Command '['jupyter', 'nbconvert', '--to', 'pdf', '--no-input', '--output', 'chronic_ACH_exposition_metrics', 'notebooks\\chronic_ACH_exposition_metrics.ipynb']' returned non-zero exit status 4294967295.

Failed to export notebook: Command '['jupyter', 'nbconvert', '--to', 'pdf', '--no-input', '--output', 'chronic_ACH_exposition_metrics', 'notebooks\\chronic_ACH_exposition_metrics.ipynb']' returned non-zero exit status 4294967295.